



GLUTACK[®]

ATRAUMATIC HERNIA MESH FIXATION SYSTEM



PIONEERING POSSIBILITY.

PRODUCT DESCRIPTION:

Glutack is a user-friendly device for the laparoscopic fixation of hernia mesh. Providing precise and consistent delivery of Glubran2 with every pull of the trigger, to minimise surgical complications and the potential pain associated with the use of tacks and staples.

ADVANTAGES:

CLINICAL BENEFITS 🕀



- · No clips or tacks: no tissue penetration, no nerve entrapment, no vessel damage
- · Less post-op pain
- Reduced post-op complications
- · Adhesive delivery controlled at all times no drips
- Tip designed to avoid clogging and sticking
- · Fixation even in high risk anatomical locations: around the Triangle of Doom, Triangle of Pain and close to the Diaphragm

STRONG (⊕)



· High shear strength equivalent to the current fixing methods: more than 9N/cm²

FAST ©



- · Very quick device preparation
- Rapid, controlled Glutack delivery
- · Each precision Glutack adheres the mesh to the tissue immediately

PRECISE Q



- Accurately controlled, repeatable Glutack volume (0.0125ml)
- No product wastage

VERSATILE (19)



- · Articulating tip system to reach difficult areas
- · Special Tip design and super rigid shaft to facilitate accurate mesh positioning

FLEXIBLE (2)



- · Multiple angles of approach facilitated by articulating tip system (0 to 90 degrees)
- · Glutacks can be delivered at any angle

INTUITIVE (4)



- Simple "point and shoot" design
- · Handle geometry provides user comfort and multiple grip
- · Audible and visual indicators confirm Glutack delivery and number of Glutacks left

GLUBRAN[®]2

Synthetic biodegradable cyanoacrylate glue, modified by the addition of a monomer synthesised by the manufacturer

Polymerises quickly in contact with live tissue and wet environment creating a thin and elastic film having high tensile strength which enables strong fixation to the tissues.

Ready for use with high adhesive, sealant haemostatic and bacteriostatic properties.

Efficient antiseptic barrier against the most diffused infective or pathogenic agents during the surgical intervention

Class III medical device authorised to be used in open, laparoscopic surgery and endovascular uses.

PRODUCT NAME **PRODUCT CODE**

GEMS-GB-DS-50

Glutack Device 50 Drops - Single - Sterile







